



Elektro

G E N E S I S

High Quality Innovations



Elektro Genesis

No 1/24 Stebondale Street, London, United Kingdom

+44 7552016275








+97 1567797033

info@elektrogenesis.com

www.elektrogenesis.com

Distributed by: Elektro Genesis ME LLC
D9 District 2 Jumeirah park Dubai UAE PO BOX 691403
STAND NO. Z4.J17

McPherson (Str.)		EG 1201
L= 9.0 cm		
Tip 0.5 mm		
McPherson (crv.)		EG 1202
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 5		EG 1203
L= 11.5 cm		
Tip 0.0 mm		
Adson		EG 1204
L= 12.1 cm		
Tip 0.1 mm		
Jeweler # 4		EG 1205
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 7		EG 1206
L= 11.5 cm		
Tip 0.5 mm		
Coaptation		EG 1207
L= 10.8 cm		
Tip 0.5 mm		

Coaptation Nadler		EG 1208
L= 10.2 cm		
Tip 0.5 mm		
Iris (Str.)		EG 1209
L= 10.2 cm		
Tip 0.5 mm		
Iris (crv.)		EG 1210
L= 10.8 cm		
Tip 1.0 mm		
Straight		EG 1211
L= 16.0 cm		
Tip 1.0 mm		
Straight		EG 1212
L= 18.0 cm		
Tip 1.0 mm		
Straight		EG 1213
L= 20.0 cm		
Tip 1.0 mm		
Straight		EG 1214
L= 22.0 cm		
Tip 1.0 mm		

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



EG 1215

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



EG 1216

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



EG 1217

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



EG 1218

Bayonet

L= 16.0 cm

Tip 1.0 mm



EG 1219

Bayonet

L= 18.0 cm

Tip 1.0 mm



EG 1220

Bayonet

L= 20.0 cm

Tip 1.0 mm



EG 1221

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



EG 1222

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



EG 1223

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



EG 1224

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



EG 1225

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



EG 1226

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



EG 1227





McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



EG 1228

McPherson (crv.) L= 11.5 cm Tip 0.5 mm		EG 1229
Jeweler # 5 L= 11.5 cm Tip 0.5 mm		EG 1230
Jeweler # 5 L= 11.5 cm Tip 0.5 mm		EG 1231
Jeweler # 7 L= 11.5 cm Tip 0.5 mm		EG 1232
Adson L= 12.1 cm Tip 1.0 mm		EG 1233
Straight L= 16.0 cm Tip 1.0 mm		EG 1234
Straight L= 18.0 cm Tip 1.0 mm		EG 1235

European Bipolar Forceps



Straight

L= 20.0 cm

Tip 1.0 mm



EG 1236

Straight

L= 22.0 cm

Tip 1.0 mm



EG 1237

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



EG 1238

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



EG 1239

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



EG 1240

Bayonet

L= 16.0 cm

Tip 1.0 mm



EG 1241








Bayonet

L= 18.0 cm

Tip 1.0 mm



EG 1242

<p>Bayonet</p> <p>L= 22.0 cm</p> <p>Tip 1.0 mm</p>		EG 1243
<p>Bayonet (Ang.Up)</p> <p>L= 16.0 cm</p> <p>Tip 1.0 mm</p>		EG 1244
<p>Bayonet (Ang.Up)</p> <p>L= 18.0 cm</p> <p>Tip 1.0 mm</p>		EG 1245
<p>Bayonet (Ang.Up)</p> <p>L= 20.0 cm</p> <p>Tip 1.0 mm</p>		EG 1246
<p>Bayonet (Ang.Down)</p> <p>L= 16.0 cm</p> <p>Tip 1.0 mm</p>		EG 1247
<p>Bayonet (Ang.Up)</p> <p>L= 18.0 cm</p> <p>Tip 1.0 mm</p>		EG 1248
<p>Bayonet (Ang.Up)</p> <p>L= 20.0 cm</p> <p>Tip 1.0 mm</p>		EG 1249

Straight

L= 16.0 cm

Tip 1.0 mm



EG 1250

Straight

L= 20.0 cm

Tip 1.0 mm



EG 1251

Bayonet

L= 16.0 cm

Tip 1.0 mm



EG 1252







Bayonet

L= 20.0 cm

Tip 1.0 mm



EG 1253

McPherson (Str.)		EG 1254
L= 9.0 cm		
Tip 0.5 mm		
McPherson (crv.)		EG 1255
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 5		EG 1256
L= 11.5 cm		
Tip 0.5 mm		
Adson		EG 1257
L= 12.1 cm		
Tip 1.0 mm		
Jeweler # 4		EG 1258
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 7		EG 1259
L= 11.5 cm		
Tip 0.5 mm		
Coaptation		EG 1260
L= 10.8 cm		
Tip 0.5 mm		

Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



EG 1261

Iris (Str.)

L= 10.8 cm

Tip 0.5 mm



EG 1262

Iris (crv.)

L= 10.8 cm

Tip 1.0 mm



EG 1263

Straight

L= 18.0 cm

Tip 1.0 mm



EG 1264

Straight

L= 20.0 cm

Tip 1.0 mm



EG 1265

Straight

L= 20.0 cm

Tip 1.0 mm



EG 1266

Straight

L= 22.0 cm

Tip 1.0 mm



EG 1267

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



EG 1268

Straight (Ang.)

L= 18.5 cm

Tip 1.0 mm



EG 1269

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



EG 1270

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



EG 1271

Bayonet

L= 16.0 cm

Tip 1.0 mm



EG 1272

Bayonet

L= 18.0 cm

Tip 1.0 mm



EG 1273

Bayonet

L= 20.0 cm

Tip 1.0 mm



EG 1274

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



EG 1275

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



EG 1276

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



EG 1277

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



EG 1278

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



EG 1279

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



EG 1280

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



EG 1281

McPherson (Crv.)
L= 9.0 cm
Tip 0.5 mm



EG 1282

Coaptation Nadler
L= 10.8 cm
Tip 0.5 mm



EG 1283

Coaptation
L= 10.8 cm
Tip 0.5 mm



EG 1284

Jeweler # 5
L= 11.0 cm
Tip 0.5 mm



EG 1285

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



EG 1286

Tenzel
L= 11.5 cm
Tip 0.5 mm



EG 1287

Jeweler # 7		EG 1289
L= 11.5 cm		
Tip 0.5 mm		
Iris (Str.)		EG 1290
L= 10.8 cm		
Tip 0.5 mm		
Iris (Crv.)		EG 1291
L= 10.8 cm		
Tip 1.0 mm		
Adson		EG 1292
L= 12.1 cm		
Tip 1.0 mm		
Semkin		EG 1293
L= 14.0 cm		
Tip 1.0 mm		
Cushing		EG 1294
L= 16.0 cm		
Tip 1.0 mm		
Cushing		EG 1295
L= 11.5 cm		
Tip 1.0 mm		

Cushing (Gursh)

L= 18.0 cm

Tip 1.0 mm



EG 1295

Semkin

Irrigation



EG 1296

Gerald Bayonet

Irrigation



EG 1297

Jansen Bayonet

Irrigation



EG 1298

Scollive Green-Wood

Irrigation



EG 1299

Hardy Bayonet

Irrigation



EG 12100

Bayonet Forceps
Irrigation



EG 12101

Bayonet Forceps
Irrigation



EG 12102

Straight Forceps
Irrigation










EG 12103

Straight Forceps
Irrigation



EG 12104

<p>McPherson (Str.)</p> <p>L= 9.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Jansen Bayonet L= 20.0 cm Tip size : 1.0 mm</p> <p>08-030-10FNS</p> <p>EG 12105</p>
<p>McPherson (Crv.)</p> <p>L= 9.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Jansen Bayonet L= 18.0 cm Tip size : 1.0 mm</p> <p>08-028-10FNS</p> <p>EG 12106</p>
<p>Coaptation Nadler</p> <p>L= 10.8 cm</p> <p>Tip 0.5 mm</p>	 <p>Cushing (Str.) L= 16.0 cm Tip size : 1.0 mm</p> <p>08-016-10FNS</p> <p>EG 12107</p>
<p>Coaptation</p> <p>L= 10.8 cm</p> <p>Tip 0.5 mm</p>	 <p>Adson L= 12.1 cm Tip size : 1.0 mm</p> <p>08-011-10FNS</p> <p>EG 12108</p>
<p>Jeweler # 5</p> <p>L= 11.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Jeweler # 5 L= 11.5 cm Tip size : 0.5 mm</p> <p>08-005-05FNS</p> <p>EG 12109</p>
<p>Jeweler # 4</p> <p>L= 11.5 cm</p> <p>Tip 0.5 mm</p>	 <p>McPherson (Crv.) L= 9.0 cm Tip size : 0.5 mm</p> <p>08-002-05FNS</p> <p>EG 12110</p>
<p>Tenzel</p> <p>L= 11.5 cm</p> <p>Tip 0.5 mm</p>	 <p>McPherson (Str.) L= 9.0 cm Tip size : 0.5 mm</p> <p>08-001-05FNS</p> <p>EG 12111</p>

Bayonet Full Non-Stick
L= 24.0 cm
Tip size : 1.0 mm

10-044-10FNS

Jeweler # 7
L= 11.5 cm
Tip 0.5 mm



EG 12112
08-008-05

Iris (Str.)
L= 10.8 cm
Tip 0.5 mm



Bayonet Full Non-Stick
L= 22.0 cm
Tip size : 1.0 mm

10-042-10FNS
EG 12113

Iris (Crv.)
L= 10.8 cm
Tip 1.0 mm



Bayonet Full Non-Stick (Irrigation)
L= 22.0 cm
Tip size : 1.0 mm

10-042-10FNS IR EG 12114

Adson
L= 12.1 cm
Tip 1.0 mm



Bayonet Full Non-Stick (Irrigation)
L= 20.0 cm
Tip size : 1.0 mm

10-040-10FNS IR EG 12115

Semkin
L= 14.0 cm
Tip 1.0 mm



Bayonet Full Non-Stick
L= 22.0 cm
Tip size : 1.0 mm

EG 12116

Cushing
L= 16.0 cm
Tip 1.0 mm



Bayonet Full Non-Stick
L= 20.0 cm
Tip size : 1.0 mm

30-040-10FNS EG 12117

Cushing
L= 11.5 cm
Tip 1.0 mm



Bayonet Full Non-Stick
L= 18.0 cm
Tip size : 1.0 mm

30-038-10FNS
EG 12118

<p>McPherson (Str.)</p> <p>L= 9.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Bayonet Full Non-Stick L= 16.0 cm Tip size : 1.0 mm</p> <p>EG 12119</p>
	30-036-10FNS
<p>McPherson (Crv.)</p> <p>L= 9.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Bayonet Full Non-Stick L= 22.0 cm Tip size : 1.5 mm</p> <p>EG 12120</p>
	10-042-15FNS
<p>Coaptation Nadler</p> <p>L= 10.8 cm</p> <p>Tip 0.5 mm</p>	 <p>Bayonet Full Non-Stick L= 20.0 cm Tip size : 1.5 mm</p> <p>EG 12121</p>
	10-040-15FNS
<p>Coaptation</p> <p>L= 10.8 cm</p> <p>Tip 0.5 mm</p>	 <p>Bayonet Full Non-Stick L= 18.0 cm Tip size : 1.0 mm</p> <p>EG 12122</p>
	10-038-10FNS
<p>Jeweler # 5</p> <p>L= 11.0 cm</p> <p>Tip 0.5 mm</p>	 <p>Straight (Ang.) Full Non-Stick L= 16.0 cm Tip size : 1.0 mm</p> <p>EG 12123</p>
	10-026-10FNS
<p>Jeweler # 4</p> <p>L= 11.5 cm</p> <p>Tip 0.5 mm</p>	 <p>Jeweler # 5 Full Non-Stick L= 11.5 cm Tip size : 0.5 mm</p> <p>EG 12124</p>
	10-005-05FNS
<p>Tenzel</p> <p>L= 11.5 cm</p> <p>Tip 0.5 mm</p>	 <p>Straight Full Non-Stick L= 12.0 cm Tip size : 1.0 mm</p> <p>EG 12125</p>
	10-012-10FNS
<p>Tenzel</p> <p>L= 11.5 cm</p> <p>Tip 0.5 mm</p>	 <p>Jeweler # 5 Full Non-Stick L= 11.5 cm Tip size : 0.5mm</p> <p>EG 12126</p>
	10-006-05FNS

Gerald Straight
Single Use



Jansen Bayonet
Irrigation



Jansen Bayonet
Single Use



Jeweler # 5
Single Use



Adson Forceps
Single Use



Cushing Forceps
Size: 16.0 cm



Jansen Bayonet ang-up
Size: 20.0 cm



Jansen Bayonet

Size: 18.0 cm



EG 12134

Jeweler # 4

Size: 11.5 cm



tip size 0.5 mm

08-004S

EG 12135

Adson Forceps

Size: 12.1 cm



tip size 1.0 mm

18-011S

EG 12136

Jansen Bayonet

Size: 18.0 cm



tip size 1.0 mm

18-028_s EG 12137

Jansen Bayonet

Size: 20.0 cm



EG 12138

Single Use European Bipolar Forceps



Bayonet Forceps
Size: 18.0 cm



EG 12139

Bayonet Forceps
Size: 20.0 cm



EG 12140

Straging Forceps
Size: 20.0 cm



EG 12141

Straight Forceps (ang)
Size: 14.0 cm



EG 12142

Straight Forceps
Size: 14.0 cm



EG 12143

Jeweler # 5
Single Use



EG 12144

Straight Forceps
Single Use



EG 12145

Bayonet Forceps
Single Use



EG 12146

Straight Forceps
Single Use



EG 12147

Straight Forceps (ang)
Single Use



EG 12148

Monopolar Cables



Silicone Monopolar Cables 10" (3.0m)

EG 12149

Instrument Connector: Screw Type

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Cables 10" (3.0m)

EG 12150

Instrument Connector: 4.0mm Male Fitting

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Generators



Silicone Monopolar Endoscopic Cable 10" (3.0m)

EG 12151

Instrument Connector: 4.0mm Female Plug

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Endoscopic Cable 10" (3.0m)

EG 12152

Instrument Connector: 4.0mm Female Plug

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Generators



Silicone Diathermy Cable 10" (3.0m)

EG 12153

Instrument Connector: 4.8mm Male Fitting

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie



Silicone Diathermy Cable 10" (3.0m)

EG 12154

Instrument Connector: 4.8mm Male Fitting

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Bipolar Cables 10" (3.0m)
EG 12155**

Instrument Connector: USA 2-Pin Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus and Mbst Bipolar Generator



**Silicone Bipolar Cables 10" (3.0m)
EG 12156**

Instrument Connector: European Flat Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10" (3.0m)
EG 12157**

Instrument Connector: European Flat Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus & Mbst Bipolar Generator



**Silicone Bipolar Cables 10" (3.0m)
EG 12158**

Instrument Connector: USA2-Pin Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10" (3.0m)
EG 12159**

Instrument Connector: USA2-Pin Plug
Generator Connector: Erbe Type
Compatibility: Erbe, Siemens, Downs



**Silicone Bipolar Cables 10" (3.0m)
EG 12160**

Instrument Connector: USA2-Pin Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley

**Silicone Bipolar Cables 10" (3.0m)
EG 12161**

Instrument Connector: European Flat Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley





**Turnur Warwick Scissor
(Straight)**
EG 12162



**Turnur Warwick Scissor
(Curved)**
EG 12163



**Surgical Scissor
(Straight)**
EG 12164



**Surgical Scissor
(Curved)**
EG 12165



**Abbey Sub Needle
17.5cm**
EG 12166



**Willson Hey Forceps
(Straight)**
EG 12167

**Richie Forceps
(Straight)**
EG 12168











**Willson Hey Forceps
(Curved)**
EG 12169



**Artery Forceps
(Straight)**
EG 12170

**Artery Forceps
(Curved)**
EG 12171



Adson		EG 12172
Mcldone		EG 12173
Gillie		EG 12174
Dissecting Forceps		EG 12175
Gerald (1x2)		EG 12176
Adson Bayonet		EG 12177
Adson Bayonet Angle Up		EG 12178
Remington Hob		EG 12179



Artery Forceps Bi-Clamp

EG 12180



Artery Forceps Non Stick

EG 12181



Artery Forceps

EG 12182



Artery Forceps

EG 12183

Silicone Bi-Clamp Connecting cord

EG 12184





Straight

EG 12185

Angled

EG 12186

Angled

EG 12187

Angled 45°

EG 12188

Angled 90°

EG 12189



Bi-Turbinate Probe

Size: 17.5cm

EG 12190

Electrodes 2.4mm



Blade Electrode
Disposable: EG-12191



Ball Electrode (3.0 mm)
Disposable: EG-12192



Circle Loop Electrode 4.0 mm
Reusable: EG-12193



Ball Electrode (crv.) (3.0 mm)
Disposable: EG-12194



Blade Electrode (crv.)
Disposable: EG-12195



Circle Loop Electrode 4.0 mm
Reusable: EG-12196



Needle Electrode
Disposable: EG-12197



Ball Electrode (4.0 mm)
Disposable: EG-12198



Circle Loop Electrode 8.0 mm
Reusable: EG-12199



Needle Electrode (crv.)
Disposable: EG-12201



Ball Electrode (crv.) (4.0 mm)
Disposable: EG-12202



Diamond Electrode 8.0 mm
Reusable: EG-12200



Ball Electrode (2.0 mm)
Disposable: EG-12203



Ball Electrode (5.0 mm)
Disposable: EG-12204



Ball Electrode (crv.) (2.0 mm)
Disposable: EG-12205



Ball Electrode (crv.) (5.0 mm)
Disposable: EG-12206



Tip Cleaner EG-12207

Holster EG-12210



Bladeelectrode (Shaft 13.0 cm)
Disposable: EG-12208



Bladeelectrode (Shaft 13.0 cm)
Disposable: EG-12209



Bladeelectrode (Shaft 13.0 cm)
Disposable: EG-12211



Bladeelectrode (Shaft 13.0 cm)
Disposable: EG-12212



Conizationelectrod
e Disposable :
EG-12213



Conization Electrode
Disposable EG-12214



Conization Electrode
Disposable EG-12215



Conization Electrode
Disposable EG-12216



Conization Electrode
Disposable EG-12217



Conization Electrode
Disposable EG-12218



Extender Electrode
EG-12219 Size: 7.0 cm



Extender Electrode
EG-12220 Size: 13.0 cm



Extender Electrode
EG-12221 Size: 18.0 cm



Blade Electrode
Reusable EG-12222



Ball Electrode
Reusable EG-12223



Blade Electrode (Crv)
Reusable EG-12224



Ball Electrode (Crv) 3mm
Reusable EG-12225



Needle Electrode
Reusable EG-12226



Ball Electrode (4mm)
Reusable EG-12227



Round Loop Electrode
(7.0cm) 4.0mm
Reusable EG-12228



Needle Electrode (Crv)
Reusable EG-12229



Ball Electrode (Crv)
(4mm) Reusable EG-12230



Round Loop Electrode
(7.0cm) 6.0mm
Reusable EG-12231



Ball Electrode (2mm)
Reusable EG-12232



Ball Electrode (5mm)
Reusable EG-12233



Round Loop Electrode
(7.0cm) 8.0mm
Reusable EG-12234



Ball Electrode (Crv) (2mm)
Reusable EG-12235



Ball Electrode (Crv) (5mm)
Reusable EG-12236



Round Loop Electrode
(7.0cm) 10.0mm
Reusable EG-12237



Blade Electrode (Shaft 13.0cm)
Reusable EG-12238



Plate Electrode Small
Reusable EG-12239



Needle Electrode (Shaft 13.0cm)
Reusable EG-12240



Ball Electrode 3mm (Shaft 13.0 cm)
Reusable EG-12241



Plate Electrode Large
Reusable EG-12242



Ball Electrode 4mm (Shaft 13.0cm)
Reusable EG-12243



Ball Electrode 5mm (Shaft 13.0cm)
Reusable EG-12244

Patient Return Plate Reusable
5m Cable EG-12245



Patient Return Plate Reusable
(5m Cable) EG-12246



Patient Return Cable
EG-12247-a (5m cable)
EG-12248-b (3m Cable)

Patient Return Cable
EG-12249 -a (5m cable)
EG-12250-b (3m Cable)



S.S Patient Return Plate
(16x12cm) EG-12251

WITHOUT SUCTION



EG-12252 L-Hook Electrode (36cm)



EG-12254 J-Hook Electrode (36cm)



EG-12256 Needle-Hook Electrode (36cm)



EG-12258 Spatula-Hook Electrode (36cm)



EG-12260 Ball-Hook Electrode (36cm)

WITH SUCTION



EG-12253 L-Hook Electrode (Irrigation) (36cm)



EG-12255 J-Hook Electrode (Irrigation) (36cm)



EG-12257 Needle-Hook Electrode (Irrigation) (36cm)



EG-12259 Spatula-Hook Electrode (Irrigation) (36cm)



EG-12261 Ball-Hook Electrode (Irrigation) (36cm)

Laparoscopic Electrodes





Graves Vaginal Speculum
with smoke evacuation tube
Small: EG-12262-a
Medium: EG-12262-b
Large: EG-12262-c



VU-Max Speculum
with smoke evacuation tube
Small: EG-12263-a
Medium: EG-12263-b
Large: EG-12263-c



Pederson Vaginal Speculum
with smoke evacuation tube
Small: EG-12264-a
Medium: EG-12264-b
Large: EG-12264-c



VU-More (Graves) Speculum
with smoke evacuation tube
Small: EG-12266-a
Medium: EG-12266-b
Large: EG-12266-c

Cusco Speculum
with smoke evacuation tube
Small: EG-12265-a
Medium: EG-12265-b
Large: EG-12265-c



Collin Speculum
with smoke evacuation tube
Small: EG-12267-a
Medium: EG-12267-b
Large: EG-12267-c





**Weisman Speculum
with Suction (Left Open)**
Small: EG-12268-a
Medium: EG-12268-b
Large: EG-12268-c



**Weisman Speculum
with Suction (Right Open)**
Small: EG-12269-a
Medium: EG-12269-b
Large: EG-12269-c



Dressing Forceps
EG-12270-a (21.0 cm)
EG-12270-b(25.0 cm)



Tissue Forceps
EG-12271-a (21.0 cm)
EG-12271-b (25.0 cm)



Iris Hook
EG-12272 (25.0 cm)



Emmet Hook
EG-12273 (25.0 cm)



Straight Hook
EG-12274 (25.0 cm)



Two-Prong Hook
EG-12275 (25.0 cm)



Three-Prong Hook
EG-12276 (25.0 cm)



Lateral Vaginal
Retractor - EG-12277



Kogan Endocervical
Speculum - EG-12278



Lateral Vaginal Retractor
(wide-Open) EG-12279



Rigby Retractor
EG-12280



Rigby Retractor
EG-12281



3mm EG-12282



4mm EG-12283



3mm EG-12284



**Kogan Endocervical
Speculum**
EG-12287



**Lateral Vaginal
Retractor**
EG-12285



**Graves Speculum
with smoking evacuation tube**
Small: EG-12286-a
Medium: EG-12286-b
Large: EG-12286-c



**Pederson Graves Speculum
with smoking evacuation tube**
Small: EG-12289-a
Medium: EG-12289-b
Large: EG-12289-c



**Collin Speculum
with smoking evacuation tube**
Small: EG-12288-a
Medium: EG-12288-b
Large: EG-12288-c





EG Headquarters

No 1/24 Stebondale Street, London, United Kingdom

+44 7552016275

+97 1567797033

info@elektrogenesis.com

www.elektrogenesis.com

Distributed by: Electro Genesis ME LLC

D9 District 2 Jumeirah park Dubai UAE PO BOX 691403

STAND NO. Z4.J17

info@elekt.com